



N Antibody, Biotin conjugated

Product Code	CSB-RA33255D1GMY
Abbreviation	N
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0DTC9
Immunogen	Recombinant Human Novel Coronavirus Nucleoprotein (N) (1-419aa) (CSB-EP3325GMY)
Species Reactivity	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	Mouse scFv fusion with human IgG1 Fc
Clonality	Monoclonal
Product Type	Recombinant Antibody
Immunogen Species	Human Novel Coronavirus (SARS-CoV-2/ 2019-nCoV)
Research Area	Microbiology
Accession NO.	1A6
Description	Recombinant anti-SARS-CoV-2 Nucleoprotein Mouse ScFv is expressed from 293 cells (HEK293) with a human IgG1 Fc tag on C-terminal.