




TNF Antibody

Product Code	CSB-PA023955KA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P01375
Immunogen	C term -peptide of Human SFPQ
Raised In	Rabbit
Species Reactivity	Human,Mouse,Rat
Tested Applications	ELISA,WB,IHC;WB:1:1000-1:2000,IHC:1:50-1:100
Form	Liquid
Storage Buffer	Store at -20oC or -80oC. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Affinity purification
Isotype	IgG
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Gene Names	TNF

Image



Western blot analysis of extracts of HeLa and NIH3T3 cells, using TNF α antibody.

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.