





CEACAM8 Antibody

Product Code	CSB-PA005168ESR2HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31997
Immunogen	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 8 protein (35-320AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Gene Names	CEACAM8

Image



Immunohistochemistry of paraffin-embedded human colon tissue using CSB-PA005168ESR2HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human lung tissue using CSB-PA005168ESR2HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.