



KDM6A Antibody, FITC conjugated

Product Code	CSB-PA025774LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15550
Immunogen	Recombinant Human Lysine-specific demethylase 6A protein (620-749AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Gene Names	KDM6A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.