




TMPRSS2 Antibody

Product Code	CSB-PA971239
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15393
Immunogen	Fusion protein of human TMPRSS2
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; ELISA:1:1000-1:2000, WB:1:200-1:1000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience?Cancer; Cell biology
Gene Names	TMPRSS2

Image



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane: A549 cells, Primary antibody: CSB-PA971239(TMPRSS2 Antibody) at dilution 1/350, Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution, Exposure time: 2 minutes